

---

# McCoy's 5A (With HEPES) Product manual

## Basic Information

Cat.NO	Size	Shelf	Form	Storage	Transportation
CMB0077	500mL	12 months	Liquid	Store at 2-8°C away from light	Room Temperature

## Product Introduction

McCoy's 5A medium was designed by Thomas McCoy and others in 1959. It was originally used specifically for the culture of Novikoff Hepatoma cells. McCoy's 5A medium is a modified Basal Medium 5A. Unlike other culture media, McCoy's 5A contains reduced glutathione, bacterial peptone, and a high concentration of glucose. McCoy's 5A is widely used for the culture of various types of primary cells, such as bone marrow, skin, gums, kidneys, spleen, lungs, rat embryos, omentum, etc. In addition, McCoy's 5A medium is also used for tissue biopsy culture, cell line establishment, and the culture of some lymphocytes and cells that are difficult to culture.

HEPES is an excellent biological buffer with no toxic effect on cells. The culture medium with HEPES added can maintain a constant pH range for a long time, which can effectively prevent the adverse effects of large pH fluctuations in the culture medium on cell growth.

## Instructions

1. Balance the culture medium and related solutions in a water bath or at room temperature, and prepare the culture medium required for the experimental cells;
2. Cell inoculation: Remove the cells to be cultured from the original culture container, wash with appropriate culture medium or PBS, and adherent cells need to be digested with trypsin;
3. Collect the cells by centrifugation, centrifuge at 1000rpm for 3 min at room temperature, and discard the supernatant;
4. Add fresh culture medium to resuspend the cells. Then inoculate the cell suspension into the culture bottle with the corresponding volume of culture medium, mix gently, and culture at 37°C and 5% CO<sub>2</sub> saturated humidity. Observe and replace fresh culture medium regularly according to cell growth and cell density.

## Precautions

1. During the entire process, be sure to pay attention to aseptic operation to avoid contamination;
2. To maintain the best use effect of this product, do not perform freeze-thaw treatment;

---

3. This product is only used for research or further research, not for diagnosis and treatment.