

Product Name: MYH1 (9C7) Mouse Monoclonal Antibody
Catalog #: AMM00722

Summary

Production Name	MYH1 (9C7) Mouse Monoclonal Antibody
Description	Primary antibody
Host	Mouse
Application	ICC/IF
Reactivity	Human,Mouse,Rat,Fruitfly,Nematode

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	MYH1
Alternative Names	fast skeletal muscle; MYH3; HEMHC; Myosin-1; Myosin heavy chain 2x; MyHC-2x; Myosin heavy chain IIx/d; MyHC-IIx/d; Myosin heavy chain; skeletal muscle; adult 1
Gene ID	4619
SwissProt ID	P12882

Application

Dilution Ratio	IF: 1/50-1/200
Molecular Weight	-

Product Name: MYH1 (9C7) Mouse Monoclonal Antibody
Catalog #: AMM00722


Background

Muscle contraction.

Research Area

Signal Transduction

Image Data



Immunofluorescence analysis of MYH1 (9C7) in Immunofluorescence Staining of nematode tissue using MYH1 (9C7) antibody.

Note

For research use only.