

**Product Name: STAT1 (9G8) Mouse Monoclonal Antibody**  
**Catalog #: AMM00787**

---

## Summary

<b>Production Name</b>	STAT1 (9G8) Mouse Monoclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Mouse
<b>Application</b>	IHC-P
<b>Reactivity</b>	Human,Rat,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	STAT1
<b>Alternative Names</b>	STAT1; Signal transducer and activator of transcription 1-alpha/beta; Transcription factor ISGF-3 components p91/p84
<b>Gene ID</b>	6772
<b>SwissProt ID</b>	P42224

## Application

<b>Dilution Ratio</b>	IHC: 1/50-1/100
<b>Molecular Weight</b>	-

**Product Name: STAT1 (9G8) Mouse Monoclonal Antibody**  
**Catalog #: AMM00787**

---


## Background

The protein encoded by this gene is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators.


## Research Area

Signal Transduction

## Image Data



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using STAT1 (9G8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using STAT1 (9G8) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## Note

For research use only.