

**Product Name: TNF alpha (6H2) Mouse Monoclonal Antibody**  
**Catalog #: AMM00887**



---

## Summary

<b>Production Name</b>	TNF alpha (6H2) Mouse Monoclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Mouse
<b>Application</b>	WB,IHC-P
<b>Reactivity</b>	Human,Rat,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	TNF
<b>Alternative Names</b>	TNF; TNFA; TNFSF2; Tumor necrosis factor; Cachectin; TNF-alpha; Tumor necrosis factor ligand superfamily member 2; TNF-a
<b>Gene ID</b>	7124
<b>SwissProt ID</b>	P01375

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100
<b>Molecular Weight</b>	Calculated MW: 26 kDa; Observed MW: 17,25 kDa

**Product Name: TNF alpha (6H2) Mouse Monoclonal Antibody**  
**Catalog #: AMM00887**




## Background

This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR2. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.


## Research Area

Immunology

## Image Data



Western blot analysis of TNF alpha in Recombinant Human TNF a Protein using TNF  $\alpha$  antibody.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil Tissue using TNF  $\alpha$  antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## Note

For research use only.

**Product Name: TNF alpha (6H2) Mouse Monoclonal  
Antibody**  
**Catalog #: AMM00887**

---

