

Product Name: SSB (4B10) Mouse Monoclonal Antibody
Catalog #: AMM03555



Summary

Production Name	SSB (4B10) Mouse Monoclonal Antibody
Description	Primary antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	SSB
Alternative Names	La; LARP3; La/SSB
Gene ID	6741
SwissProt ID	P05455

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 50 kDa

Background

Binds to the 3' poly(U) terminii of nascent RNA polymerase III transcripts, protecting them from exonuclease digestion and

Product Name: SSB (4B10) Mouse Monoclonal Antibody
Catalog #: AMM03555




facilitating their folding and maturation.

Research Area

Immunology

Image Data



Western blot analysis of La Antigen in MCF-7, A549, Jurkat and HeLa lysates using La Antigen antibody.

Note

For research use only.