

Product Name: GRP78 BiP (9C7) Mouse Monoclonal Antibody
Catalog #: AMM03619



Summary

Production Name	GRP78 BiP (9C7) Mouse Monoclonal Antibody
Description	Primary antibody
Host	Mouse
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	HSPA5 HSPA5; GRP78; 78 kDa glucose-regulated protein; GRP-78; Endoplasmic reticulum
Alternative Names	luminal Ca(2+)-binding protein grp78; Heat shock 70 kDa protein 5; Immunoglobulin heavy chain-binding protein; BiP
Gene ID	3309
SwissProt ID	P11021

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 72 kDa; Observed MW: 78 kDa

Product Name: GRP78 BiP (9C7) Mouse Monoclonal Antibody
Catalog #: AMM03619




Background

When Chinese hamster K12 cells are starved of glucose, the synthesis of several proteins, called glucose-regulated proteins (GRPs), is markedly increased. Hendershot et al. (1994) (PubMed 8020977) pointed out that one of these, GRP78 (HSPA5), also referred to as 'immunoglobulin heavy chain-binding protein' (BiP), is a member of the heat-shock protein-70 (HSP70) family and is involved in the folding and assembly of proteins in the endoplasmic reticulum (ER).


Research Area

Tags & Cell Markers

Image Data




Western blot analysis of GRP78 BiP (9C7) in HeLa, rat Liver lysates using GRP78/Bip antibody.




Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma using GRP78/Bip antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.


Product Name: GRP78 BiP (9C7) Mouse Monoclonal Antibody
Catalog #: AMM03619



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using GRP78/Bip antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Western blot analysis of GRP78 BiP in Jurkat, MCF-7, HeLa lysates using GRP78 BiP (9C7) antibody



Immunocytochemistry analysis of GRP78 BiP in HeLa using Bip/GRP78 (9C7) antibody.

Product Name: GRP78 BiP (9C7) Mouse Monoclonal Antibody
Catalog #: AMM03619



Immunoprecipitation analysis of GRP78 BiP in HeLa lysates using Bip/GRP78 (9C7) antibody.

Note

For research use only.