

---

## Summary

<b>Production Name</b>	Importin 9 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ICC/IF,IP
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	IPO9
<b>Alternative Names</b>	importin 9; Imp9; Ran-binding protein 9; RanBP9
<b>Gene ID</b>	55705
<b>SwissProt ID</b>	Q96P70

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
<b>Molecular Weight</b>	Calculated MW: 116 kDa; Observed MW: 116 kDa

## Background

---

**Product Name: Importin 9 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01433**




Nuclear transport receptor that mediates nuclear import of proteins, such as histones, proteasome and actin.


## Research Area

Signal Transduction

## Image Data



Immunocytochemistry analysis of Importin 9 (green) in A549 using Importin 9 antibody, and DAPI (blue).



Western blot analysis of Importin 9/RANBP9 in Jurkat, HeLa lysates using Importin 9/RANBP9 antibody.

## Note

For research use only.