

Product Name: PAH Rabbit Monoclonal Antibody
Catalog #: AMRe01458



Summary

Production Name	PAH Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PAH
Alternative Names	Phe-4-monooxygenase
Gene ID	5053
SwissProt ID	P00439

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Background

Product Name: PAH Rabbit Monoclonal Antibody
Catalog #: AMRe01458




Catalyzes the hydroxylation of L-phenylalanine to L-tyrosine.

Research Area

Signal Transduction

Image Data



Western blot analysis of PAH in 3T3, C6, rat Brain, K562 lysates using PAH antibody.

Note

For research use only.