

Summary

Production Name	SIRT2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	SIRT2
Alternative Names	SIRT2; SIR2L; SIR2L2; NAD-dependent protein deacetylase sirtuin-2; Regulatory protein SIR2 homolog 2; SIR2-like protein 2
Gene ID	22933
SwissProt ID	Q8IXJ6

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43 kDa

Product Name: SIRT2 Rabbit Monoclonal Antibody
Catalog #: AMRe01521




Background

Participate in a variety of cellular functions, including histone deacetylation, gene silencing, chromosomal stability, and aging. SIRT2, a human homolog of the yeast SIR2, functions as transcriptional silencing mediator at mating-type loci, telomeres and ribosomal gene clusters. SIRT2 expression increases dramatically during mitosis and is multiply phosphorylated at the G(2)/M transition of the cell cycle.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of SIRT2 in K562, HeLa lysates using SIRT2 antibody.

Note

For research use only.