

Product Name: R Cadherin Rabbit Monoclonal Antibody
Catalog #: AMRe01528



Summary

Production Name	R Cadherin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CDH4
Alternative Names	cadherin 4; type 1; R-cadherin (retinal); CAD4; RCAD; R-CAD
Gene ID	1002
SwissProt ID	P55283

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 100 kDa; Observed MW: 130 kDa

Background

Product Name: R Cadherin Rabbit Monoclonal Antibody
Catalog #: AMRe01528




Cadherins are calcium-dependent cell adhesion proteins.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of R-Cadherin/Cadherin 4 in mouse brain, U87-MG lysates using R-Cadherin antibody.



Western blot analysis of R-Cadherin/Cadherin 4 in rat Brain, C6 lysates using R-Cadherin/Cadherin 4 antibody.

Note

For research use only.