
Summary

Production Name	USP11 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	USP11
Alternative Names	UHX1; USP11
Gene ID	8237
SwissProt ID	P51784

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 110 kDa; Observed MW: 110 kDa

Background

Product Name: USP11 Rabbit Monoclonal Antibody
Catalog #: AMRe01529




Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Inhibits the degradation of target proteins by the proteasome. Plays a role in the regulation of pathways leading to NF-kappa-B activation. Plays a role in the regulation of DNA repair after double-stranded DNA breaks.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of USP11 in Jurkat, rat Brain, HeLa lysates using USP11 antibody.

Note

For research use only.