

Product Name: ALDH1L1 Rabbit Monoclonal Antibody
Catalog #: AMRe01624



Summary

Production Name	ALDH1L1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ALDH1L1
Alternative Names	FDH; FTHFD; 10-fTHF; 10-FTHFDH
Gene ID	10840
SwissProt ID	O75891

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 99 kDa; Observed MW: 99 kDa

Background

Product Name: ALDH1L1 Rabbit Monoclonal Antibody
Catalog #: AMRe01624




Highly expressed in liver, pancreas and kidney.

Research Area

Neuroscience

Image Data



Western blot analysis of ALDH1L1 in rat Brain, 3T3 lysates using ALDH1L1 antibody.

Note

For research use only.