

Product Name: ARHGAP25 Rabbit Monoclonal Antibody
Catalog #: AMRe01674



Summary

Production Name	ARHGAP25 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ARHGAP25
Alternative Names	Rho GTPase activating protein 25; KAIA0053
Gene ID	9938
SwissProt ID	P42331

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Background

Product Name: ARHGAP25 Rabbit Monoclonal Antibody
Catalog #: AMRe01674




GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Research Area

Signal Transduction

Image Data



Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.

Note

For research use only.