

Summary

Production Name	BAG3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	BAG3
Alternative Names	BAG3; BIS; BAG family molecular chaperone regulator 3; BAG-3; Bcl-2-associated athanogene 3; Bcl-2-binding protein Bis; Docking protein CAIR-1
Gene ID	9531
SwissProt ID	O95817

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 80 kDa

Product Name: BAG3 Rabbit Monoclonal Antibody
Catalog #: AMRe01711




Background

Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Has anti-apoptotic activity.

Research Area

Cell Biology

Image Data



Western blot analysis of Bag3 in K562, rat Brain lysates using BAG3 antibody.

Note

For research use only.