

Product Name: CD55 Rabbit Monoclonal Antibody
Catalog #: AMRe01797



Summary

Production Name	CD55 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CD55
Alternative Names	CD55; CR; DAF; Complement decay-accelerating factor; CD55
Gene ID	1604
SwissProt ID	P08174

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 78 kDa

Background

Product Name: CD55 Rabbit Monoclonal Antibody
Catalog #: AMRe01797




This protein recognizes C4b and C3b fragments that condense with cell-surface hydroxyl or amino groups when nascent C4b and C3b are locally generated during C4 and c3 activation. Interaction of daf with cell-associated C4b and C3b polypeptides interferes with their ability to catalyze the conversion of C2 and factor B to enzymatically active C2a and Bb and thereby prevents the formation of C4b2a and C3bBb, the amplification convertases of the complement cascade.

Research Area

Immunology

Image Data



Western blot analysis of CD55 in HeLa lysates using CD55 antibody.

Note

For research use only.