

Product Name: Chromogranin A Rabbit Monoclonal Antibody
Catalog #: AMRe01825



Summary

Production Name	Chromogranin A Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CHGA
Alternative Names	beta Granin; Chromogranin A; ChrA; SP-I; Vasostatin I; Vasostatin-2; Vasostatin II; Pancreastatin; CHGA
Gene ID	1113
SwissProt ID	P10645

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 92 kDa

Product Name: Chromogranin A Rabbit Monoclonal Antibody
Catalog #: AMRe01825




Background

Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

Research Area

Neuroscience

Image Data



Western blot analysis of Chromogranin A in PC-12 lysates using Chromogranin A antibody.

Note

For research use only.