

Product Name: Collagen VI alpha 1 Rabbit Monoclonal Antibody
Catalog #: AMRe01841



Summary

Production Name	Collagen VI alpha 1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	COL6A1
Alternative Names	OPLL; BTHLM1; UCHMD1
Gene ID	1291
SwissProt ID	P12109

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Molecular Weight	Calculated MW: 109 kDa; Observed MW: 147 kDa

Product Name: Collagen VI alpha 1 Rabbit Monoclonal Antibody
Catalog #: AMRe01841




Background

Collagen VI acts as a cell-binding protein.

Research Area

Cell Biology

Image Data



Western blot analysis of Collagen VI in Human placenta lysates using Collagen VI alpha 1 antibody.

Note

For research use only.