
Summary

Production Name	CT A Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PCYT1A
Alternative Names	CCT A; CCT alph; CTPCT; PCYT1; Pcyt1a
Gene ID	5130
SwissProt ID	P49585

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa

Background

Product Name: CT A Rabbit Monoclonal Antibody
Catalog #: AMRe01856




Controls phosphatidylcholine synthesis.

Research Area

Neuroscience

Image Data



Western blot analysis of CT A in MCF-7, HeLa lysates using CT A antibody.

Note

For research use only.