
Summary

Production Name	Follistatin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FST
Alternative Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene
Gene ID	10468
SwissProt ID	P19883

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Background

Product Name: Follistatin Rabbit Monoclonal Antibody
Catalog #: AMRe01995




Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Research Area

Signal Transduction

Image Data



Western blot analysis of Follistatin in rat Brain lysates using Follistatin antibody.

Note

For research use only.