

**Product Name: Glutathione Reductase Rabbit
Monoclonal Antibody
Catalog #: AMRe02042**

Summary

Production Name	Glutathione Reductase Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	GSR
Alternative Names	GR; HEL-75; HEL-S-122m
Gene ID	2936
SwissProt ID	P00390

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 51 kDa

**Product Name: Glutathione Reductase Rabbit
Monoclonal Antibody**
Catalog #: AMRe02042


Background

Maintains high levels of reduced glutathione in the cytosol.

Research Area

Signal Transduction

Image Data



Western blot analysis of Glutathione Reductase in 3T3, HeLa lysates using Glutathione Reductase antibody.

Note

For research use only.