
Summary

Production Name	HBEGF Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	HBEGF
Alternative Names	DTR; DTS; DTSF; HB-EGF; HBEGF; HEGFL
Gene ID	1839
SwissProt ID	Q99075

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Product Name: HBEGF Rabbit Monoclonal Antibody
Catalog #: AMRe02071




May be involved in macrophage-mediated cellular proliferation. It is mitogenic for fibroblasts and smooth muscle but not endothelial cells.

Research Area

Signal Transduction

Image Data



Western blot analysis of HBEGF in K562, C6, 3T3, HeLa lysates using HBEGF antibody.

Note

For research use only.