

Summary

Production Name	hnRNP K Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	HNRNPK
Alternative Names	HNRNPK; HNRPK; Heterogeneous nuclear ribonucleoprotein K; hnRNP K; Transformation up-regulated nuclear protein; TUNP
Gene ID	3190
SwissProt ID	P61978

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 58 kDa

Product Name: hnRNP K Rabbit Monoclonal Antibody
Catalog #: AMRe02107




Background

Facilitate pre-mRNA processing and transport of mRNA from the nucleus to cytoplasm. hnRNP K contains three unique structural motifs termed KH domains that bind poly(C) DNA and RNA sequences. Intricate architecture enables hnRNP K to facilitate mRNA biosynthesis, transcriptional regulation, and signal transduction. Research studies have shown that cytoplasmic hnRNP K expression is increased in oral squamous cell carcinoma and pancreatic cancer, and may be a potential prognostic factor.


Research Area

Epigenetics and Nuclear Signaling

Image Data




Western blot analysis of hnRNP K in 3T3, MEF lysates using hnRNP K antibody.




Western blot analysis of hnRNP K in C6, CHO-K1, HeLa lysates using hnRNP K antibody.

Product Name: hnRNP K Rabbit Monoclonal Antibody
Catalog #: AMRe02107



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using hnRNP K antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of hnRNP K (green) in MCF-7 using hnRNP K antibody, and DAPI (blue).

Note

For research use only.