

**Product Name: Interferon gamma Receptor 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02136**



Summary

Production Name	Interferon gamma Receptor 1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	IFNGR1
Alternative Names	IFNGR1; Interferon gamma receptor 1; IFN-gamma receptor 1; IFN-gamma-R1; CDw119; CD antigen CD119
Gene ID	3459
SwissProt ID	P15260

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 100 kDa

**Product Name: Interferon gamma Receptor 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02136**




Background

Receptor for interferon gamma. Two receptors bind one interferon gamma dimer.

Research Area

Immunology

Image Data



Western blot analysis of IFNGR1 in K562, C6, 3T3, HeLa lysates using Interferon gamma Receptor 1 antibody.

Note

For research use only.