

Product Name: Integrin beta 5 Rabbit Monoclonal Antibody
Catalog #: AMRe02165

Summary

Production Name	Integrin beta 5 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Hamster,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ITGB5
Alternative Names	Integrin beta-5
Gene ID	3693
SwissProt ID	P18084

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 88 kDa

Product Name: Integrin beta 5 Rabbit Monoclonal Antibody
Catalog #: AMRe02165


Background

Integrin alpha-V/beta-5 (ITGAV:ITGB5) is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand. (Microbial infection) Integrin ITGAV:ITGB5 acts as a receptor for adenovirus type C.

Research Area

Signal Transduction

Image Data



Western blot analysis of Integrin beta 5 in C6, CHO-K1, HeLa lysates using Integrin beta 5 antibody.

Note

For research use only.