

**Product Name: Phospho-LAT (Tyr200) Rabbit
Monoclonal Antibody
Catalog #: AMRe02207**

Summary

Production Name	Phospho-LAT (Tyr200) Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	LAT
Alternative Names	LAT; Linker for activation of T-cells family member 1; 36 kDa phospho-tyrosine adapter protein; pp36; p36-38
Gene ID	27040
SwissProt ID	O43561

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 28 kDa; Observed MW: 36-38 kDa

**Product Name: Phospho-LAT (Tyr200) Rabbit
Monoclonal Antibody
Catalog #: AMRe02207**


Background

Required for TCR (T-cell antigen receptor)- and pre-TCR-mediated signaling, both in mature T-cells and during their development. Involved in FCGR3 (low affinity immunoglobulin gamma Fc region receptor III)-mediated signaling in natural killer cells and FCER1 (high affinity immunoglobulin epsilon receptor)-mediated signaling in mast cells. Couples activation of these receptors and their associated kinases with distal intracellular events such as mobilization of intracellular calcium stores, PKC activation, MAPK activation or cytoskeletal reorganization through the recruitment of PLCG1, GRB2, GRAP2, and other signaling molecules.

Research Area

Immunology

Image Data



Western blot analysis of Phospho-LAT (Tyr200) in Jurkat lysates using Phospho-LAT (Tyr200) antibody.

Note

For research use only.