
Summary

Production Name	LGALS14 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	LGALS14
Alternative Names	Charcot-Leyden crystal protein 2; CLC2; Galectin-14; Gal-14
Gene ID	56891
SwissProt ID	Q8TCE9

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Product Name: LGALS14 Rabbit Monoclonal Antibody
Catalog #: AMRe02212




Binds beta-galactoside and lactose.


Research Area

Cell Biology

Image Data



Immunocytochemistry analysis of LGALS14 (green) in K562 using LGALS14 antibody, and DAPI (blue).



Western blot analysis of LGALS14 in C6 lysates using LGALS14 antibody.

Note

For research use only.