

Product Name: MMP26 Rabbit Monoclonal Antibody
Catalog #: AMRe02265



Summary

Production Name	MMP26 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Hamster,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MMP26
Alternative Names	Matrix metalloproteinase-26; Endometase; Matrilysin-2
Gene ID	56547.0
SwissProt ID	Q9NRE1

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 19 kDa

Background

Product Name: MMP26 Rabbit Monoclonal Antibody
Catalog #: AMRe02265




May hydrolyze collagen type IV, fibronectin, fibrinogen, beta-casein, type I gelatin and alpha-1 proteinase inhibitor.

Research Area

Signal Transduction

Image Data



Western blot analysis of MMP26 in Jurkat, C6, CHO-K1 lysates using MMP26 antibody.

Note

For research use only.