
Summary

Production Name	NDUFS6 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	NDUFS6
Alternative Names	CI-13kA; CI13KDA; CI-13kD-A
Gene ID	4726
SwissProt ID	O75380

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 14 kDa; Observed MW: 14 kDa

Background

Product Name: NDUFS6 Rabbit Monoclonal Antibody
Catalog #: AMRe02323




Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of NDUFS6 in rat heart lysates using NDUFS6 antibody.



Western blot analysis of NDUFS6 in COS7, Raji, 3T3, Jurkat lysates using NDUFS6 antibody.

Note

For research use only.