
Summary

Production Name	SLC12A1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	SLC12A1
Alternative Names	BSC1; NKCC2
Gene ID	6557
SwissProt ID	Q13621

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 121 kDa; Observed MW: 150 kDa

Background

Product Name: SLC12A1 Rabbit Monoclonal Antibody
Catalog #: AMRe02346




Electrically silent transporter system. Mediates sodium and chloride reabsorption. Plays a vital role in the regulation of ionic balance and cell volume.

Research Area

Signal Transduction

Image Data



Western blot analysis of SLC12A1/NKCC2 in HeLa lysates using SLC12A1 antibody.

Note

For research use only.