
Summary

Production Name	Nova1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	NOVA1
Alternative Names	Nova-1
Gene ID	4857
SwissProt ID	P51513

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Background

Product Name: Nova1 Rabbit Monoclonal Antibody
Catalog #: AMRe02352




May regulate RNA splicing or metabolism in a specific subset of developing neurons.

Research Area

Neuroscience

Image Data



Western blot analysis of Nova1 in rat lung lysates using Nova1 antibody.

Note

For research use only.