
Summary

Production Name	NUDT9 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	NUDT9
Alternative Names	NUDT10
Gene ID	53343
SwissProt ID	Q9BW91

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 39 kDa

Background

Product Name: NUDT9 Rabbit Monoclonal Antibody
Catalog #: AMRe02365




Hydrolyzes ADP-ribose (ADPR) to AMP and ribose 5'-phosphate.


Research Area

Signal Transduction

Image Data



Western blot analysis of NUDT9 in K562, HeLa lysates using NUDT9 antibody.



Immunocytochemistry analysis of NUDT9(green) in HeLa using NUDT9 antibody, and DAPI(blue)

Note

For research use only.