
Summary

Production Name	PEX11B Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PEX11B
Alternative Names	PEX14B; PEX11-BETA
Gene ID	8799
SwissProt ID	O96011

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 28 kDa; Observed MW: 28 kDa

Background

Product Name: PEX11B Rabbit Monoclonal Antibody
Catalog #: AMRe02427




Involved in peroxisomal proliferation (PubMed:9792670). May regulate peroxisome division by recruiting the dynamin-related GTPase DNM1L to the peroxisomal membrane (PubMed:12618434). Promotes membrane protrusion and elongation on the peroxisomal surface (PubMed:20826455).


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of PEX11B in T47D, HL-60 lysates using PEX11B antibody.



Western blot analysis of PEX11B in 3T3 lysates using PEX11B antibody.

Note

For research use only.