

Product Name: PKC gamma Rabbit Monoclonal Antibody
Catalog #: AMRe02450

Summary

Production Name	PKC gamma Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PRKCG
Alternative Names	PRKCG; PKCG; Protein kinase C gamma type; PKC-gamma
Gene ID	5582
SwissProt ID	P05129

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 78 kDa; Observed MW: 78 kDa

Product Name: PKC gamma Rabbit Monoclonal Antibody
Catalog #: AMRe02450


Background

This is a calcium-activated, phospholipid-dependent, serine- and threonine-specific enzyme.


Research Area

Signal Transduction

Image Data



Western blot analysis of PKC gamma in rat brain lysates using PKC gamma antibody.



Western blot analysis of PKC in HeLa lysates using PKC antibody.

Note

For research use only.