

Product Name: Prealbumin Rabbit Monoclonal Antibody
Catalog #: AMRe02470



Summary

Production Name	Prealbumin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TTR
Alternative Names	Prealbumin; TBPA; Amyloidosis I; Transthyretin
Gene ID	7276
SwissProt ID	P02766

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Product Name: Prealbumin Rabbit Monoclonal Antibody
Catalog #: AMRe02470




Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Research Area

Neuroscience

Image Data



Western blot analysis of Prealbumin in rat Brain lysates using Prealbumin antibody.

Note

For research use only.