
Summary

Production Name	STAT6 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	STAT6
Alternative Names	STAT6; Signal transducer and activator of transcription 6; IL-4 Stat
Gene ID	6778
SwissProt ID	P42226

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 94 kDa; Observed MW: 110 kDa

Background

Product Name: STAT6 Rabbit Monoclonal Antibody
Catalog #: AMRe02641




STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

Research Area

Signal Transduction

Image Data



Western blot analysis of STAT6 in Jurkat, CHO-K1, HeLa lysates using STAT6 antibody.

Note

For research use only.