

Product Name: PLEKHA1 Rabbit Monoclonal Antibody
Catalog #: AMRe02665



Summary

Production Name	PLEKHA1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PLEKHA1
Alternative Names	PH domain-containing family A member 1
Gene ID	59338
SwissProt ID	Q9HB21

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 46 kDa; Observed MW: 46 kDa

Background

Product Name: PLEKHA1 Rabbit Monoclonal Antibody
Catalog #: AMRe02665




Binds specifically to phosphatidylinositol 3,4-diphosphate (PtdIns3,4P2), but not to other phosphoinositides. May recruit other proteins to the plasma membrane.

Research Area

Signal Transduction

Image Data



Western blot analysis of TAPP1 in K562 lysates using PLEKHA1 antibody.

Note

For research use only.