

Product Name: T Box Protein 3 Rabbit Monoclonal Antibody
Catalog #: AMRe02671



Summary

Production Name	T Box Protein 3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TBX3
Alternative Names	T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL
Gene ID	6926
SwissProt ID	O15119

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 79 kDa; Observed MW: 79 kDa

Product Name: T Box Protein 3 Rabbit Monoclonal Antibody
Catalog #: AMRe02671


Background

Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Tbx3 in 3T3 lysates using T Box Protein 3 antibody.

Note

For research use only.