

Summary

Production Name	TCTP Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TPT1
Alternative Names	Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23
Gene ID	7178
SwissProt ID	P13693

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Product Name: TCTP Rabbit Monoclonal Antibody
Catalog #: AMRe02676




Background

Involved in calcium binding and microtubule stabilization.


Research Area

Immunology

Image Data



Western blot analysis of TCTP in K562, rat Brain, C6, 3T3, HeLa lysates using TCTP antibody.



Immunocytochemistry analysis of TCTP(green) in HeLa using TCTP antibody, and DAPI(blue)

Note

For research use only.