

**Product Name: Transmembrane Protein 192 Rabbit
Monoclonal Antibody
Catalog #: AMRe02696**

Summary

Production Name	Transmembrane Protein 192 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TMEM192
Alternative Names	Transmembrane protein 192
Gene ID	201931
SwissProt ID	Q8IY95

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 31 kDa; Observed MW: 31 kDa

**Product Name: Transmembrane Protein 192 Rabbit
Monoclonal Antibody
Catalog #: AMRe02696**


Background

Enables protein homodimerization activity. Located in several cellular components, including late endosome; lysosomal membrane; and perinuclear region of cytoplasm.

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of TMEM192 in K562, HeLa lysates using Transmembrane Protein 192 antibody.

Note

For research use only.