

**Product Name: Transmembrane Protein 85 Rabbit
Monoclonal Antibody
Catalog #: AMRe02699**

Summary

Production Name	Transmembrane Protein 85 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	EMC4
Alternative Names	PIG17; TMEM85
Gene ID	51234
SwissProt ID	Q5J8M3

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

**Product Name: Transmembrane Protein 85 Rabbit
Monoclonal Antibody**
Catalog #: AMRe02699


Background

May mediate anti-apoptotic activity.


Research Area

Cell Biology

Image Data



Western blot analysis of TMEM85 in 3T3 lysates using Transmembrane Protein 85 antibody.



Western blot analysis of TMEM85 in Jurkat, rat Brain, C6, CHO-K1, HeLa lysates using TMEM85 antibody.

Note

For research use only.