
Summary

Production Name	Tensin 4 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TNS4
Alternative Names	CTEN; C-terminal tensin-like protein
Gene ID	84951
SwissProt ID	Q8IZW8

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 77 kDa; Observed MW: 90 kDa

Background

Product Name: Tensin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe02702




May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton .
May promote apoptosis, via its cleavage by caspase-3.


Research Area

Signal Transduction

Image Data



Western blot analysis of TNS4/CTEN in Raw264.7, C6 lysates using Tensin 4 antibody.



Western blot analysis of TNS4/CTEN in K562 lysates using TNS4/CTEN antibody.

Note

For research use only.