
Summary

Production Name	TRAF6 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TRAF6
Alternative Names	TRAF6; RNF85; TNF receptor-associated factor 6; E3 ubiquitin-protein ligase TRAF6; Interleukin-1 signal transducer; RING finger protein 85
Gene ID	7189
SwissProt ID	Q9Y4K3

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 60 kDa; Observed MW: 60 kDa

Product Name: TRAF6 Rabbit Monoclonal Antibody
Catalog #: AMRe02707




Background

TRAFs (TNF receptor-associated factors) are a family of multifunctional adaptor proteins that bind to surface receptors and recruit additional proteins to form multiprotein signaling complexes capable of promoting cellular responses. Members of the TRAF family share a common carboxy-terminal TRAF domain which mediates interactions with associated proteins; many also contain amino-terminal Zinc/RING finger motifs.

Research Area

Cell Biology

Image Data



Western blot analysis of TRAF6 in 293 lysates using TRAF6 antibody.

Note

For research use only.