

**Product Name: Translin Associated Protein X Rabbit
Monoclonal Antibody
Catalog #: AMRe02712**



Summary

Production Name	Translin Associated Protein X Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TSNAX
Alternative Names	TRAX
Gene ID	7257
SwissProt ID	Q99598

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 33 kDa; Observed MW: 33 kDa

**Product Name: Translin Associated Protein X Rabbit
Monoclonal Antibody
Catalog #: AMRe02712**




Background

Acts in combination with TSN as an endonuclease involved in the activation of the RNA-induced silencing complex (RISC).
Possible role in spermatogenesis.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TRAX in 3T3, mouse testis lysates using Translin Associated Protein X antibody.



Western blot analysis of TRAX in K562, rat Brain, C6, HeLa lysates using TRAX antibody.

Note

For research use only.