

Product Name: ULK3 Rabbit Monoclonal Antibody
Catalog #: AMRe02745



Summary

Production Name	ULK3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ULK3
Alternative Names	Serine/threonine-protein kinase ULK3; Ulk3; unc 51 like kinase 3 (C. elegans); Unc-51-like kinase 3
Gene ID	25989
SwissProt ID	Q6PHR2

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Product Name: ULK3 Rabbit Monoclonal Antibody
Catalog #: AMRe02745




Background

Serine/threonine protein kinase which enhances GLI1 and GLI2 transcriptional activity and consequently positively regulates GLI-dependent SHH signaling. May exert this function by promoting GLI1 nuclear localization. Phosphorylates in vitro GLI2, as well as GLI1 and GLI3, although less efficiently.

Research Area

Signal Transduction

Image Data



Western blot analysis of ULK3 in K562, HeLa lysates using ULK3 antibody.

Note

For research use only.