

Product Name: Hex Rabbit Monoclonal Antibody
Catalog #: AMRe02806



Summary

Production Name	Hex Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	HHEX
Alternative Names	HEX; HHEX; HMPH; PRH; PRHX
Gene ID	3087
SwissProt ID	Q03014

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 30 kDa

Background

Product Name: Hex Rabbit Monoclonal Antibody
Catalog #: AMRe02806




Recognizes the DNA sequence 5'-ATTAA-3'. Transcriptional repressor. May play a role in hematopoietic differentiation. Establishes anterior identity at two levels; acts early to enhance canonical WNT-signaling by repressing expression of TLE4, and acts later to inhibit NODAL-signaling by directly targeting NODAL.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Hex in Jurkat, HeLa lysates using Hex antibody.

Note

For research use only.