
Summary

Production Name	IFRD1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	IFRD1
Alternative Names	PC4; TIS7
Gene ID	3475
SwissProt ID	O00458

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa

Background

Product Name: IFRD1 Rabbit Monoclonal Antibody
Catalog #: AMRe02812




Could play a role in regulating gene activity in the proliferative and/or differentiative pathways induced by NGF.

Research Area

Cell Biology

Image Data



Western blot analysis of IFRD1 in Ramos lysates using IFRD1 antibody.

Note

For research use only.