
Summary

Production Name	RBM43 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RBM43
Alternative Names	C2orf38
Gene ID	375287
SwissProt ID	Q6ZSC3

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Background

Product Name: RBM43 Rabbit Monoclonal Antibody
Catalog #: AMRe02817




Predicted to enable RNA binding activity.

Research Area

-

Image Data



Western blot analysis of RBM43 in HeLa lysates using RBM43 antibody.

Note

For research use only.