

Product Name: Elongation Factor Ts Rabbit Monoclonal Antibody
Catalog #: AMRe02818

Summary

Production Name	Elongation Factor Ts Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TSFM
Alternative Names	EFTS; EFTSMT
Gene ID	10102
SwissProt ID	P43897

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Product Name: Elongation Factor Ts Rabbit Monoclonal Antibody
Catalog #: AMRe02818




Background

Associates with the EF-Tu.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TFSM in HeLa, A549, HL-60, U2OS lysates using Elongation Factor Ts antibody.

Note

For research use only.