
Summary

Production Name	CD229 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	LY9
Alternative Names	hly9; mLY9; CD229; SLAMF3
Gene ID	4063
SwissProt ID	Q9HBG7

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 72 kDa; Observed MW: 72 kDa

Background

Product Name: CD229 Rabbit Monoclonal Antibody
Catalog #: AMRe02933




Self-ligand receptor of the signaling lymphocytic activation molecule (SLAM) family.

Research Area

-

Image Data



Western blot analysis of CD229 in HeLa, A549, HL-60, U251 lysates using CD229 antibody.

Note

For research use only.