

**Product Name: Elongation Factor 1A2 Rabbit
Monoclonal Antibody
Catalog #: AMRe02978**

Summary

Production Name	Elongation Factor 1A2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	EEF1A2
Alternative Names	HS1; STN; EF1A; STNL; MRD38; EEF1AL; EIEE33; EF-1-alpha-2
Gene ID	1917
SwissProt ID	Q05639

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa

**Product Name: Elongation Factor 1A2 Rabbit
Monoclonal Antibody
Catalog #: AMRe02978**


Background

This protein promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of EEF1A2 in mouse brain, rat brain, mouse heart lysates using Elongation Factor 1A2 antibody.



Western blot analysis of EEF1A2 in HeLa, A549 lysates using EEF1A2 antibody.

Note

For research use only.